

Self Seal Sterilisation Pouches



LOADING INSTRUMENT INTO POUCH



PEELING OFF STRIP FROM ADHESIVE BAND



FOLDING PAPER FLAP ALONG PERFORATION



PRESSING DOWN FIRMLY TO CREATE SEAL



PEELING OPEN FROM OPPOSITE END AFTER PROCESSING

- ▶ Manufactured by Granton Medical in Sheffield, England in a Class 10,000 Clean Room facility on unique, high speed, fully validated plant.
- ▶ Heavyweight 70gsm paper
- ▶ Blue tinted transparent 2-ply laminate film allows easy viewing of the contents of the pouch
- ▶ Suitable for both steam and ethylene oxide sterilisation
- ▶ Each pouch printed with Class A process indicators for both steam and ethylene oxide sterilisation
- ▶ Each pouch printed with a lot number for full manufacturing traceability
- ▶ Self Seal strip eliminates the need for heat seal equipment and provides a validated seal
- ▶ Off-Set webs for easy opening
- ▶ Tack seals to prevent curling and reduce external contamination
- ▶ Comprehensive size range, custom sizes available by arrangement

Standard sizes available from stock:-

CODE	NOMINAL SIZE	BOX QTY	CASE QUANTITY
11110	2¼" x 4" / 60mm x 100mm	Box 200	10 boxes of 200 per case
11120	3½" x 5½" / 90mm x 140mm	Box 200	10 boxes of 200 per case
11128	3½" x 8" / 90mm x 200mm	Box 200	10 boxes of 200 per case
11130	3½" x 9" / 90mm x 230mm	Box 200	10 boxes of 200 per case
11140	5¼" x 10" / 135mm x 255mm	Box 200	10 boxes of 200 per case
11144	5¼" x 14" / 135mm x 360mm	Box 200	10 boxes of 200 per case
11150	7½" x 12" / 190mm x 305mm	Box 200	5 boxes of 200 per case
11153	7½" x 13" / 190mm x 330mm	Box 200	5 boxes of 200 per case
11185	12" x 15" / 305mm x 381mm	Box 200	5 boxes of 200 per case

Granton Medical
Quality Healthcare Products

Granton Ragg Ltd, Parkway Close, Parkway Industrial Estate, Sheffield S9 4WJ England

Tel: +44 (0)114 275 7290 Fax: +44 (0)114 272 7853

E-mail: info@granton.co.uk Website: www.granton-medical.co.uk