

Granton® SILVER

Surgical Blades & Scalpels



- ▶ A compact range of scalpels from Granton®, Knifemakers since 1601
- ▶ Packaging designed to minimise waste and reduce cost
- ▶ Consistent quality and reliable performance for professional users
- ▶ Single use products with CE mark accreditation
- ▶ Sterile products sterilised by gamma irradiation

BLADE PATTERN	STAINLESS STERILE BLADES (100's)	CARBON STERILE BLADES (100's)	CARBON NON-STERILE BLADES (100's)	STERILE DISPOSABLE SCALPELS (10's)
10	91010	91110	91410	91310
10A	91000	91100	91400	91300
11	91011	91111	91411	91311
12	91012	91112	91412	91312
15	91015	91115	91415	91315
20	91020	91120	91420	91320
21	91021	91121	91421	91321
22	91022	91122	91422	91322
23	91023	91123	91423	91323
24	91024	91124	91424	91324


METAL SCALPEL HANDLES - Regular alloy and stainless handles accurately machined and finished to BS and ISO standards to ensure easy blade fitment and interchangeability with other brands of blades. Supplied non-sterile and also available in plastic.

HANDLE PATTERN	ALLOY (10's)	STAINLESS (10's)	PLASTIC (100's)	BLADE FITMENT
3	90203	90213	90303	Small/No.3
7	90207	90217	-	Small/No.3
9	90209	90219	-	Small/No.3
3L	90223	90233	-	Small/No.3
4	90204	90214	90304	Large/No.4
4L	90224	90234	-	Large/No.4

SURGICAL BLADES - Manufactured from hardened and tempered surgical steel to BS 2982:1992 and BS EN 27740:1992 / ISO 7740:1985 standards. Precision ground to ensure fine cutting edges and reliable fitment with handles. Supplied in sterile foil packs or non-sterile packs for single use.

DISPOSABLE SCALPELS - Single piece plastic mouldings with removable guards to protect blades before use and safeguard against accidental cross-infection after use. Supplied in sterile paper/film packs for single use.

 Granton
Medical
Quality Healthcare Products



Granton Ragg Ltd, Parkway Close, Parkway Industrial Estate, Sheffield S9 4WJ England

Tel: +44 (0)114 275 7290 Fax: +44 (0)114 272 7853

E-mail: info@granton.co.uk Website: www.granton.co.uk

Granton®



Cert. No. LRQ 0942127