

Dental Laboratory Knives & Spatulas

A long standing leader in its field, Granton[®] produces an extensive range of knives and associated products to comply with all the latest regulations.

All our knife products are made in Sheffield, to the highest standards from the finest materials and are fully guaranteed. The Granton[®] range of Dental knives and spatulas are precision made for safe and comfortable professional use. Stainless steel blades with handles insert moulded in polypropylene to comply with best hygiene and safety practise, withstanding sterilisation by autoclave or gamma irradiation.

Dental Plaster Knives

Code No.	Blade Length
49000	2½" / 60mm
49010	3" / 75mm
49020	3" / 75mm (with Flask Opener)



Dental Wax Knives

Code No.	Blade Length
49100	Large Dental Wax Knife, Buffed & Bent
49110	Small Dental Wax Knife, Buffed & Bent
49120	Large Dental Wax Knife, Dished & Polished
49130	Small Dental Wax Knife, Dished & Polished



Dental Plaster Mixers

Code No.	Blade Length
49200	4¼" / 110mm Broad Curved
49210	4¼" / 110mm Broad Flat
49220	3⅞" / 100mm Medium Curved
49230	3⅞" / 100mm Medium Flat
49240	4¼" / 110mm Narrow Curved
49250	4¼" / 110mm Narrow Flat



Dental Laboratory Knives

Code No.	Blade Length
49400	1¾" / 45mm Clip Point
49410	1¾" / 45mm Sheep Foot



Granton® Palette Knives/Spatulas (with Flexible Blades)

Granton® Palette Knives are manufactured in the traditional way; by fully grinding the entire length of the blade. The test when choosing a Palette Knife is to press the blade onto a firm surface. It should bend from the middle of the blade, not near the handle. Keep pushing the blade so it is at 90 degrees to the handle.

A Granton® Palette knife has the flexibility to allow it to return to a dead straight "rest" position time and time again.

From our 3" Artists Palette Knife up to our 12" Spatula Palette, the feedback we receive from our customers tells us that the Granton® Palette knife range is ahead of the rest! For safe professional use in Laboratory, Science and Education establishments, Granton® blades are made to BS and ISO standards from fully hardened and tempered stainless steel, taper ground and hand finished with a non-reflective sateen surface for excellent flexibility and comfortable use under strong lighting.

Code No.	Blade Length	Description
49309	3" / 75mm - Tapered Blade	"Artists" Palette
49310	4" / 100mm - Tapered Blade	"Artists" Palette
49311	6" / 150mm - Tapered Blade	"Artists" Palette
45284	4" / 100mm	Spatula Palette
45286	6" / 150mm	Spatula Palette
45288	8" / 200mm	Spatula Palette
45280	10" / 250mm	Spatula Palette
45282	12" / 300mm	Spatula Palette

